



Mobile Medical Assistant

MPAD-800

8" Rugged Medical Tablet



Features

- Verizon Certified LTE with Android 6.0
- Type-A & C USB for Easy Connection
- Water & Drop Proof
- 10 Hours Battery Life
- Slim & Light Weight Design
- 4th Edition (UL/EN60601-1)

Specifications

Main Specifications

Processor	ARM® Cortex® - A17 Quad-Core up to 1.6 GHz
System Memory	2 GB SDRAM
OS	Android 6.0.1
LCD Display	8" (1280 x 800) TFT LCD with 2-point P-cap
Communication	WiFi b/g/n
	Bluetooth 4.0 (EDR) with integrated class 1 PA and BLE Support
I/O Port	USB Type C x 1
	USB Type A x 1
	Volume Up/Down
	3.5 mm Stereo Headphone Jack x 1
	DC-In Jack x 1
	Built-In Speaker & Microphone
Storage	nanoSIM Card Slot
	Micro SD card slot
Storage	16 GB eMMC Flash
Camera	Front 2MP, Rear 8MP with Flash LED
Sensor	Built-in Light Sensor
	Built-in G-Sensor
	Built-in Vibration Motor

Mechanical


Color	White
Dimension (W x H x D)	136.6 x 219 x 21 mm (5.38 x 8.58 x 0.83")
Gross Weight	0.67 Kg (1.48 lb)

Environmental

Operating Temperature	0°C ~ 40°C
Storage Temperature	-20°C ~ 60°C
Operating Humidity	10 ~ 75%
Water/ Dust Protection Rating	IP54

Power Supply

AC/DC Adapter	AC Power Input: 100 ~ 240 V, 50 ~ 60 Hz, 36 W
	Output: 12 V, 3 A
Battery	Lithium-Polymer High Capacity Battery



Remark:

Docking function and docking station is optional for special order.

