

Non-Powered Medication Delivery Cart with Drawers

OP-TD-SV43-13X0



Features

- Innovative medication delivery cart for the point of care
- Non-powered with add-on power options available
- Three choices of display mounts
- Auto-locking drawer system
- Truly ergonomic design
- Lightweight cart with a small footprint
- Customize with storage and accessories
- Ships fully assembled

Specifications

Main Specifications

LCD Pivot	Storage compartment for larger computers can mount to rear of cart with Universal CPU Holder accessory (#80-105-064)
Compact wheel-to-wheel base dimensions	17" x 21" (43 x 53,3 cm)
Cart weight	Varies by drawer configuration
Shipping dimensions	Varies by drawer configuration
Shipping weight	Varies by drawer configuration
Products include	LCD Pivot Worksurface CPU lock Secure CPU compartment Four-port USB hub Scanner holder bracket Mouse holder Wrist rest Back-tilt keyboard tray Left/right mouse tray Rear handle Storage basket Cable storage Four casters (two locking) StyleLink SKY software 5-year warranty for mechanical and 1-year for electrical components

Non-Powered Documentation Cart with Optional Drawers


Several Storage Options and Accessories Available:



SV Storage Bin



Several drawer configurations



SV Front Tray



SV Sani-wipes Holder



SV Wire Basket, Large

Product Information:

	Typical LCD Size	Weight Capacity	Lift	Tilt	Pan	Rotation	VESA	
Non-Powered SV40	LCD Pivot	≤ 24"	35 lbs (15,9 kg)	Lift Pivot 25 (63,5 cm)	Arm 65° Lift Pivot 25°	Arm 180° Lift Pivot 24°	Lift Pivot 90° P/L	MIS-D
		Includes secure CPU Compartment						
	5-year warranty for mechanical components							



- ① Supports displays up to 24" and provides 25 degrees of tilt and 90° rotation
 - LCD Pivot: LCD and worksurface adjust 20"(51 cm) in tandem;
 - LCD Pivot adjusts 5"(13 cm) independently;
 - Maximum LCD height adjustment = 27" (68 cm)
- ② Scanner holder bracket
- ③ Spacious worksurface: 22.4"W x 19.75"D (56,7W x 50,2D cm)
- ④ Auto-locking drawer system offers simple and efficient medication management at the point-of-care; security system stores variablelength PINs for up to 1,000 users with StyleLink SKY software
- ⑤ Back-tilt keyboard tray, wrist rest and mouse holder; Four amber LEDs provide light over keyboard tray
- ⑥ Rear storage basket keeps task-related items close for easy access
- ⑦ T-slot channel interface on the column allows you to integrate drawers and accessories
- ⑧ Four 4" (10,2 cm) casters; two front casters are lockable
- ⑨ Wheel-to-wheel base dimensions: 15.5"x 17" (39,4 x 43,2 cm)

