


Medica Display for Medical Device

# MEDDP-822

22" Slim Medical Display



## Features

- 21.5" LCD Panel with LED Backlight
- High Brightness 500nits
- Compact and TerseDesign
- DVI x 1, HDMI x 1, DPx1
- True Flat Screen with Capacitive Touch
- Speaker built-in x 2
- Physical OSD key bottom

## Specifications

### DISPLAY

Size	21.5"
Resolution	1920 x 1080 pixels
Max. Colors	16.7M
Brightness (Typ.)	500 nits
Viewing Angle (CR=10)(Typ.)	178°(H)/178°(V)
Contrast Ratio (Typ.)	1000:1
Touch Screen	Projected Capacitive Multi Touch Screen

### SYSTEM

OS Support	Windows® 10
------------	-------------

### I/O

Input Signal	VGA x 1, DVI-D x 1, HDMI x 1
USB Port	USB 2.0 x 1 for Touch Screen
Speaker	Built in 5W speaker x 2
Audio	Line-in x 1
Power	DC 12V Input x 1

### MECHANICAL AND ENVIRONMENTAL

VESA	75 x 75mm & 100 x 100mm
IP Rating	Front IP65, Back cover IPX1
Power Requirement	External power adapter, AC IN: 100 ~ 240V,12V
Operating Temperature	0°C ~ 35°C(32°F ~ 95°F)
Storage Temperature	-20°C ~ 60°C(-4°F ~ 140°F)
Storage Humidity	10% to 90% RH, noncondensing
Dimension	526 x 324 x 54mm
Packing Size	640 x 470 x 210mm
Gross Weight	6.8kg
Net Weight	5.0kg

Certifications	CE: EN 60601-1-2:2015(V4.0), EN 60601-1:2006/A1:2013/A12:2014 (V3.1) FCC: Part 18 Class B UL: ANSI/AAMI ES60601-1:2012 (V3.1) cUL: CAN/CSA-C22.2 No. 60601-1:2014 (V3.1)
----------------	--

## Ordering Information

- MEDDP-822-G1-A1-0010  
Monitor. 22".500 nits.DC 12V.Glass.HDMI cable.Medical Adapter
- MEDDP-822-P1-A1-0010  
Monitor. 22".500 nits.DC 12V.P-CAP.HDMI cable.Medical Adapter