




Medica Display for Medical Device

MEDDP-722

22" Slim Medical Display



Features

- 21.5" LCD Panel with LED Backlight
- Compact and TerseDesign
- DVI x 1, HDMI x 1, DPx1
- True Flat Screen with Capacitive Touch
- Speaker -2 x built-in
- Support 2M Camera

Specifications

DISPLAY

| | |
|-----------------------------|---|
| Size | 21.5" |
| Resolution | Full HD 1920 x 1080 |
| Max. Colors | 16.7M |
| Brightness (Typ.) | 250 nits |
| Viewing Angle (CR=10)(Typ.) | 178°(H)/178°(V) |
| Contrast Ratio (Typ.) | 3000:1 |
| Touch Screen | Projective capacitive touch screen (optional) |

SYSTEM

| | |
|------------|-------------------------|
| OS Support | Windows® 7, Windows® 10 |
| Webcam | 2MP |

I/O

| | |
|-------------|------------------------------|
| Video Input | DVIx 1, HDMI x 1, VGA x1 |
| Audio | Audio-in x 1 |
| USB | USB x 1 for touch controller |
| Speaker | 5W speaker x 2 on back side |

MECHANICAL AND ENVIRONMENTAL

| | |
|-----------------------|---|
| VESA | 75/100 mm VESA mount |
| Power Requirement | DC12V |
| Power Consumption | Full Loading :28.7Watts |
| Operating Temperature | 5°C ~ 40°C(41°F ~ 104°F) |
| Storage Temperature | -10°C ~ 60°C(14°F ~ 140°F) |
| Storage Humidity | 10%~90%@35°C, non-condensing |
| IP Rating | Protection: front IP54 - rear IPX1 |
| Dimension | 546 x 351 x 56 mm |
| Packing Size | 570 x 478 x 260 mm |
| Gross Weight | 8.4kg (18.51lb) |
| Net Weight | 5.7kg (12.56lb) |
| Certifications | CE: EN 60601-1-2:2015(V4.0), EN 60601-1:2006/ A1:2013/A12:2014(V3.1) FCC: Part 18 Class B |

Optional Accessories

- 125530036A
(AOH)(TF)Medical Power Adaptor.AC/DC.100 ~ 240V.12V.3A.36W.Adapter. ATM036T-P120(Level VI)
- 1700252001
(AOH)(TF)DVI-D Cable.(24+1)P DVI (M).(24+1)P DVI (M).200cm
- 1704191800
(AOH)(TF)HDMI Cable.HDMI 19P(F) Molding.HDMI 19P(M) Molding.180cm
- 1704151800
(AOH)(TF)VGA Cable.D-Sub 15P(M) Molding.D-Sub 15P(M) Molding.180cm
- 1700031802
AOH)(TF)Aidio Cable.3.5mm 3Pole Stereo Plug(Black).3.5mm 3Pole Stereo Plug(Green).180cm

Ordering Information

- MEDDP-722-G1-A1-0010
Monitor. 22".250 nits.DC 12V.Glass
- MEDDP-722-G1-A2-0010
Monitor. 22".250 nits.DC 12V.Glass.2M Camera
- MEDDP-722-P1-A1-0010
Monitor. 22".250 nits.DC 12V.PCT
- MEDDP-722-P2-A1-0010
Monitor. 22".250 nits.DC 12V.PCT.2M Camera