


Medical Cart Computer

Venus-244P

24" Dual Hot Swappable Battery Medical AI Cart Computer



R003-220254
T D220163003



Features

- Intel Arrow Lake Core Ultra 5 (15th Gen)
- Medical AI Cart Computer
- Increase security (TPM/Onyx screensaver)
- Imprivata compliant RFID reader (optional)
- Windows 11 and IGEL Compliant
- WiFi 7
- Dual Hot Swappable Batteries
- Easy cleanable and anti-bacterial
- Remotely Manageable by central IT departments (ORION)

Specifications

MAIN SPECIFICATIONS

Processor	Intel® Arrow Lake Core™ Ultra 5 225U
System Memory	DDR5 up to 48GB
OS Support	Microsoft® Windows 11 Ubuntu 24.04 LTS IGEL (thin client solution)
Storage	m.2 SSD up to 2TB
Wireless Communication	802.11 a/b/g/n/ac/ax (optional) WiFi 6E + Bluetooth 5.3 (optional) Wifi 7 + Bluetooth 5.4 (optional)
Touch	PCT /AG/Privacy / Filter (optional)
Speaker	3W x2
Trusted Platform Module	TPM 2.0
Security	Smart Card Reader (optional),Imprivata RFID reader (optional),Barcode reader (optional)
Waterproof	Front IP65 ; rear IPX1

DISPLAY

Display Size	23.8"
Resolution	FHD 1920 x 1080
Max. Colors	16.7M
Contrast Ratio	3000:1
Luminance (cd/m2)	250

I/O

USB	Rear I/O :USB 3.0 x2, USB Type C x1 Front button I/O : USB 3.0 x2
Video Out	HDMI 1.4 x1
COM	COM port x2 (3Kv isolation optional)
Audio	Combo audio port x1
Ethernet	1.5kv isolated Gigabit LAN x2
DC-in	19V DC-in x1

MECHANICAL AND ENVIRONMENTAL

VESA	75/100 mm VESA mount
Operating Temperature	0°C ~ 35°C(32°F ~ 95°F)
Storage Temperature	-20°C ~ 60°C(-4°F ~ 140°F)
Storage Humidity	10%~95%@40°C, non-condensing
Dimension	560(L) x 78.4(W) x 373(H)
Package Size	730(L) x 240(W) x 550(H) mm
Gross Weight	11kg(24.3lbs)
Net Weight	7kg
Certifications	CE: EN 60601-1-2: 2015 + A1: 2021 (V4.1), EN 60601-1: 2007 + A1:2013 +A2:2021 (V3.2) EN 55032:2015/A1:2020 EN55035:2017/A11:2020(ITE) IEC62368-1:2020+A11:2020 (ITE) FCC : Part 15B/ Part 18 UL: ANSI AAMI ES60601-1:2005/A1:2012/A2:2021(V3.2) cUL: CAN/CSA-C22.2 No. 60601-1:14 (IEC 60601-1: 2005+A1:2012+A2:2020,MOD) (V3.2) TELEC JAPAN TELEC

ALL PRODUCT, PRODUCT SPECIFICATIONS AND DATA ARE SUBJECT TO CHANGE WITHOUT NOTICE TO IMPROVE RELIABILITY, FUNCTION OR DESIGN OR OTHERWISE.